




# Aerosol SPA system

Pneumatic Purging-Mist Mode

## Instruction Use For Preven-Co Nasal Irrigator



### Components:

- a. Tubing
- b. Pressure Plug
- c. Irrigation Cover
- d. Irrigation Cover O-Ring Seal
- e. Pressure Connector
- f. Orifice for Mist Operating Valve
- g. Saline Solution Tank
- h. Saline solution Siphon Tubing
- i. Nasal Adapter/Mist Nozzle Holder
- j. Mist Nozzle
- k. Orifice of Mist Nozzle
- l. Tubing Adapter

### How To Install & Use: (Refer to Left Figure)



1. Unscrew the irrigation cover on saline solution tank to open.
2. Pour saline solution into the saline solution tank (20ml Max)
3. Tightly screw the cover/nasal irrigator onto saline solution tank and connect with tubing
4. Connect the other end of tubing with tubing adapter's pressure connector.
5. Connect tubing adapter's pressure plug to filtration cartridge's pressure air outlet.
6. Connect filtration cartridge's pressure plug to vacuum pump/Compressor.
7. Complete the installation of nasal cavity irrigator with vacuum pump/compressor.
8. Switch on the power and position the mist nozzle near (0-5mm) to nostril, press the orifice for mist operating valve by single finger to start the rinsing work.

## **Cleaning And Maintenance:**

- 1. Let the mist nozzle keep spray ten sec. or more fully evacuate the solution and avoid the nozzle being blocked by remaining crystal salt.**
- 2. Switch off the power then remove the filtration cartridge and tubing adapter from vacuum pump/compressor and disconnect two ends of tubing.**
- 3. Discard all unused saline solution then to clean all components (except tubing and filtration cartridge) with purified water and put them at cool place 1-2 hours for drying.**
- 4. For further sterilization, to steam all components (except tubing and filtration cartridge) with high temperature or put the components into boiled water boiling for ten minutes.**
- 5. Tubing and filtration cartridge do not need to sterilize, but, to replace the tubing every 3-6 months and filtration cartridge every two months is a suggestion.**

## **Cautions:**

- 1. The nasal irrigator is applicable for all aged users. But, for baby or younger children, the operation is suggested with the help of adult.**
- 2. The nasal irrigator is designed with non-pressure and constant spray, which is safely applicable while nasal congestion. User can swallow and spit out the solution during cleansing process.**
- 3. Besides special solution or normal saline solution, other salty water is not allowed to avoid the block on mist nozzle.**
- 4. The saline solution for contact lens, which contains chemical materials, will cause uncomfortable feeling on nasal cavity.**
- 5. One set of components for one user only. Do not share the components with other user to prevent transmitting flu. virus infection,**
- 6. For keeping best effect on operation, to replace the components every 3-6 months is a suggestion.**
- 7. Do not self-detach the irrigation cover or mist nozzle. Such action will cause block or damage.**
- 8. When spray is seriously blocked or unsuccessful, please contact the store where you brought it for repair. Do not use any tool for self-repair.**